Neuroimaging advances in Parkinson's Disease; The application of Ultra-High Field MRI

Stijn Michielse PhD
Department of Neurosurgery
Maastricht University
@stijnimaging








UTAP

- Understanding
 - 9.4T post-mortem
- Tracking
 - TRACK-PD, 7T
- Adjustment
 - sensors and adaptive DBS
- Parkinson's Disorder







Background UTAP

- ➤ Early diagnosis of Parkinson's Disease (PD)
 - Challenging / often not immediately recognised (1)
 - Heterogeneous disorder (2, 3)
 - The underlying aetiology is poorly understood (4)


1. Hughes et al. (2002); 2. Berg et al. (2014); 3. Lewis et al. (2005); 4. Thenganatt et al. (2014)







Early biomarker; swallow-tail disappears



Dorsolateral Substantia Nigra pars compacta







Work Package 1 – Understand PD

- Focus on changes in the microcircuit of the basal ganglia and brainstem
- Diffusion Weighted Imaging
- White matter visualisation and quantification
- Quantitative T2 mapping







Brain container - development














Brain container









Brain samples

- ➤ Post-mortem PD
- > Fluorinert
- ➤ Brains from UK biobank
 - > 2-3 years old, fixed in formaldehyde PBS solution









MRI protocol

- **❖**T2
 - 0.25mm isotropic GRadient Echo; entire hemisphere
 - Six echo's
 - 4 6.98ms, 11ms, 16.21ms, 20.23ms 24.46ms and 30ms
- Diffusion Weighted Imaging
 - 1mm isotropic; entire hemisphere
 - 48 random directions b-value 5009s/mm²
 - ❖ 5 low b-value volumes at 279s/mm²
 - ❖ Field of view; 144x132x180mm, TR=450ms







Shimming


- Field inhomogeneities (B0)
- ❖ Material introduced → disturbances
- Higher field strength=more inhomogeneities
- Affects image quality
- ❖B1 shim for optimal tissue flip angle
- Complex procedure at 9.4T
- Offline kT-points 3th order shim







Field inhomogeneities T2 (short echo)









Scanning


- T2 weighted; 3h (2x)
- Low b-value diffusion; 20min (5x)
- 2h per set of four directions (12x)
- Total 32h (including scan cooldowns)
- Raw data reconstruction (500Gb of data)
- Berkeley Advanced Reconstruction Toolbox (BART) in MatLab








T2 weighted – 250µm resolution







T2 weighted – 250μm resolution










Quantitative T2* map









Diffusion Weighted Imaging


- Create contrast, based on b-values (gradient strength)
- Goal; white matter quantification
- ❖White matter tracking → tractography
- Focus on limbic system and basal ganglia





White matter bundles











Terms used

- Diffusion weighted imaging
 - DWI → the MRI acquisition
- Diffusion tensor imaging
 - DTI → tensors, the smallest element
- Tractography
 - Connecting the tensors; create fibers
- Anisotropic diffusion
 - Diffusion equal to all sides







Diffusion Weighted Imaging


Noisy data

Directional dependence

Strong gradients

Multiple directions

■ B0 image → bright








DTI – color coding

Standard DTI color coding

X = red; medio-lateral

Y = green; anterior-posterior

Z = blue; inferior-superior










Reconstruct tensor

- >Smallest unit that can be measured
- ➤ Directional dependence











Isotropic movement





vs Anisotropic










2D image of tractography










Fiber tracking










Fiber tracking (2)











Pre and post central gyrus



Motor cortex



Sensory cortex







Use of the brain shell - dedicated coil









Thanks for your attention

Credits to:

Yasin Temel

Ali Jahanshahianyar

Alard Roebroeck

Jackson Boonstra

Mark Kuijf

Albert Leentjens

Amée Wolters

Margot Heijmans

Ben Poser

Heidi Jacobs

Collaborators:

Nikos Priovoulos

Koen Schruers

Liesbet Goossens

Samantha Baldi

Arash Aghamohammadi Sereshki

Stanislau Hrybouski

Wojciech Pietrasik

Nikolai Malykhin

Peter Seres

Graham Murray

Iris Lange

Jindra Bakker

Wolfgang Viechbauer

Machteld Marcelis

Jim van Os

Therese van Amelsvoort















